


Modules on Disk 9

- So we are nearly able to recover the modules operation with the links in, except that we are also seeing the following:



- Discontinuities are seen on modules on the Disk.
- It is occurring at $\sim 0.24fC$
- It is independent of trigger rate, optical power or whether the Redundancy Links are in or not.

Noise Occ. Scan with no Redundancy Links

- Compared with the SystemTest results where there is little evidence of any discontinuity.
- Threshold scan is also more uniform.

